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organized in a special body so that a user working with a UV-Cabinet is protected against UV light. Recirculator increases the maximum density of UV light making it sufficiently effective for DNA/RNA inactivation. The UV-recirculator processes 100 UV-Cabinet volumes per hour, creating permanent aseptic conditions of operation inside the UV-Cabinet.

Specially assigned moving tables T-4 (with wheel locks) with a drawer are available on request.

Advantages of Biosan UV-Cabinets:

Ozone free high density UV decontamination

Long living UV lamps (8,000 hours average)

Automatic switch off of UV-lamps when the protective screen is opened

Bactericidal flow-type recirculator providing permanent decontamination inside UV -cabinet during operation

Shockproof glass walls

Low noise, low energy consumption

Tables for installation of UV-Cabinets

UV-Cabinets with the bactericidal UV cleaner-recirculator AR is the patented Biosan solution


Specification

Walls material: Plexiglas: Polymethyl methacrylate ALTUGLAS EX

Working surface material: Steel with chemicals resistant powder coating

Open UV-lamp: 1 x 25W built in bactericidal, TUV25WG13 UV-C

Radiation type: UV ($\lambda = 253.7$ nm),

ozone-free

UV radiation level: 15 mW / cm² / sec

Digital time setting of direct UV exposure: 0–24 hrs / non-stop (increment 1 min)

UV-recirculator: 1x25W (efficiency >99% per 1 hour)

Daylight lamp (for working area illumination): 1 x TLD-15W

Thickness of side panels: 4 mm

Thickness of upper front panel: 8 mm

Thickness of screen: 8 mm

Optical transmission: 92%

UV protection: >99.90% Polymethyl methacrylate ALTUGLAS EX

Working area: 650 × 475 mm

Safety: Automatic open UV-lamp switch off when screen is open

Power outlets inside the unit: Inlet for power cords

Nominal operating voltage: 100-240 V, 50/60Hz

Power consumption: 67 W

Overall dimensions: 690 × 535 × 555 mm

Optional table: T-4 (800 × 600 × 750 mm)

Weight (net / gross): T-4 (800 × 600 × 750 mm)

Weight (net / gross): 23 / 33 kg

2/ MSC-6000, Centrifuge/Vortex Multispin

Centrifuge/vortex Multi-Spin MSC-6000 is product of many years evolution of Spin-Mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing (mix) and collecting the reagents again from the walls and cap of

the microtube (second spin). We named this repetitive algorithm of operation that is aimed at reducing the mistakes during sample preparation for PCR analysis a “sms–algorithm”.

Multi–Spin is a fully automatic device for reproducing sms–algorithm for 12 tubes at one time, thus saving time considerably. A must-have instrument for PCR and DNA analyses laboratory.

Multi Spin is four devices combined in one:

Centrifuge — Maximum RCF: up to $2350 \times g$

Vortex (3 mixing modes — soft, medium, hard; regulated time; Vortexing regulation timer 1–20 sec)

Centrifuge/Vortex;

SMS–cyclers for realization of the “sms–algorithm”.

Comparison of FVL-2400N, MSC-3000 and MSC-6000

Multi–Spin allows considerable time saving compared to Combi–Spin by automatically performing cycling program of sample mixing and spinning according to the set spin–mix–spin cycle for 12 microtubes simultaneously.

Specification

Speed control range: 1000–6000 rpm (increment 100 rpm)

RCF max: $2350 \times g$

Spin timer: 1 sec–30 min

Vortexing intensity: Soft, medium, hard

Vortexing time: 0–20 sec (increment 1 sec)

SMS-cycle regulation 1 - 999 cycles

Display: LCD, 2 x 16 signs

Safety: Lid lock

Overall dimensions:190x235x125 mm

Weight: 2.5 kg

Input current/power consumption: 24 V,
24 W (1 A)

External power supply: Input AC 100–
240 V; 50/60 Hz; Output DC 24 V

3/ FisherBrand* Electrophoresis gel
Systems, sub-gel midi-plus deal for
restriction fragment analysis, sample
prep or checking of high numbers of
samples.

- Supplied with 150mm x 70mm,
150mm x 100mm and 150mm x 150mm
gel trays
- Run up to 210 samples
- Low buffer volumes
- Multichannel pipette compatible
combs for fast gel loading
- Injection moulded

Sub-gel Midi-Plus includes all 3 tray
sizes for optimum versatility and
value. Electrophoresis, Blotting and
Immunodetection

Technical Specification - General

Dimensions [l x w], mm: 150 x 70, 150
x 100, 150 x 150 (gel)

Dimensions [l x w x h], mm: 175 x 265
x 90 (unit)

Capacity: 70 samples, max. (100 x
70mm tray)

140 samples, max. (100 x
100mm tray)

210 samples, max. (150 x
150mm tray)

Volume, mL: 500 (buffer)

No. of samples: 1, 2, 4, 10, 10MC, 12,
14MC, 16, 18MC, 20, 28MC, 30MC,
35

Thickness, mm: 0.75, 1, 1.5, 2
MC=Multichannel pipette compatible
Accessories



4/ Thermo Scientific Nalgene Labtop cooler for 32 (0.5 to 2.0 mL) tubes (holds -15°C for 2 hours)

-
- An excellent alternative to messy ice buckets for keeping reagents and enzymes cool
- Lids feature a gridded write-on surface for sample identification
- Stackable design saves valuable lab space

A convenient alternative to messy ice buckets. Body is molded from durable polycarbonate and filled with a nontoxic insulating gel; lid is made of clear polycarbonate (models 44001-10 and -25 have gel-filled white polycarbonate lid). Lids feature a gridded write-on surface for sample identification.

Models 44001-00, -10, and -25 include adapters to accommodate 0.5-mL tubes. Models 44001-20 and -30 hold tubes up to 125-mm tall when lid is latched. All models have a stainless steel handle for a secure hold during transportation.

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